

One vision, Two sharp eyes  
with Our Innovation

# CASIA2

Cornea / Anterior Segment OCT

## CASIA2 SPECIFICATIONS

### [Measurement Performance]

Resolution	Axial (Depth)	10μm or less (in tissue)
	Transverse	30μm or less (in air)
Scan range	Depth	13mm
	Transverse	Radial Scan: φ16mm
		Raster Scan: 12mm×12mm

### [Main Unit]

Scan rate	50,000 A scans / second
Stroke range of moving section	40mm(Y axis); 88mm(X axis); 45mm(Z axis)
Stroke range of chin rest	70mm
Dimensions and Weight	530(W)×560(D)×455(H)mm Approx. 33kg
Type of light source	Swept source Laser
Wavelength	1,310nm
Output power	Less than 6mW
Laser Class	Class 1

### [Power source]

Voltage	100~240V AC
Frequency	50 / 60Hz
Power consumption	170VA

### [External HDD]

Capacity	8TB or more
----------	-------------

### [Touch panel LCD monitor]

Display	Touch panel LCD monitor 20 inches or larger
---------	---

### [Workstation computer]

OS	Windows®8.1 or 10 64bit
CPU	Intel® Core i5
Memory	8GB or more
SSD	128GB
HDD	8TB or more
Data output	Printer (LAN/USB)

### Tomey Corporation [Asia-Pacific]

2-11-33 Noritakeshinmachi  
Nishi-ku, Nagoya, 451-0051, Japan  
Tel: ++81-52-581-5327  
Fax: ++81-52-561-4735  
E-mail: intl@tomey.co.jp

### Tomey GmbH [Europe]

Wiesbadener Straße 21  
90427 Nürnberg, Germany  
Tel: ++49-911-9385462-0  
Fax: ++49-911-9385462-20  
E-mail: info@tomey.de

For more information, visit our web site <http://www.tomey.com>

Always read and follow the instructions for use. Not all products, services, or offers are approved or offered in every market. Please note that the current status of approval for the labeling, instructions, and contents of the brochure may vary from one country to another.  
©2016 Tomey Corporation. Specifications are subject to change without notice. Any products mentioned herein are registered trademarks of their respective owners.

210331

One vision, Two sharp eyes with Our Innovation

# CASIA2


Cornea / Anterior Segment OCT

Efficient support of cataract surgery

- Advanced Imaging  
Deeper, wider and clearer
- Testing application for  
cataract surgery
- Fulfilling Analysis Function  
Trend analysis  
Lens shape analysis



- Advanced Imaging  
Deeper, wider and clearer
- Testing application for cataract surgery
- Fulfilling Analysis Function  
Trend analysis  
Lens shape analysis





One vision, Two sharp eyes  
with Our Innovation

# CASIA2

Cornea / Anterior Segment OCT

## Efficient support of cataract surgery


The next-generation Cornea / Anterior Segment OCT “CASIA2”, will advance cataract surgery

## Advanced Imaging


Cornea and lens shown in one image

### Deeper Imaging


With CASIA2, the light source of coherency functions is improved, and higher sensibility toward depth is realized compared to our former model. By using this new technology, it is possible to measure to a depth of 13mm from anterior cornea to posterior lens with one shot.




SS-1000



CASIA2



SS-1000




CASIA2

### Wider Imaging


Capturing images around the angle is possible in corneal topography mode. As with corneal shape analysis, it is possible to extract and analyze the angle and observe the IOL, which enables testing without switching measuring modes.

*\*For lens shape measurement, LensBiometry capture is necessary.*



### Clearer Imaging

By scanning 16 images simultaneously, clearer images are obtained.



## CASIA IOL Cataract Surgery CICS

The testing application for cataract surgery, CICS, is installed in the CASIA2, which effectively supports cataract surgery. There are two types in CICS: “Pre-op testing” and “Post-op testing”. To use their functions, capture the image within each testing protocol.

Pre-op Testing  
Pre-op  
Cataract


Screening  
IOL Cal.  
Toric IOL

Post-op Testing  
Post-op  
Cataract

Thumb  
Image


## Fulfilling analysis functions

### Trend Analysis



The color code map shows parameter changes of each corneal shape. Additionally, the simplicity of the graph means information is instinctively easy to grasp making this analysis useful for observing the keratoconus.

### Lens Shape Analysis



While capturing anterior cornea to posterior lens, it is possible to measure corneal curvature, thickness and tilt of the anterior / posterior lens.

